

Organic Application Note

Carbon, Hydrogen, and Nitrogen in Coke

Accessories 502-186 Tin Foil Cup

Calibration Standard NIST, BCR or other suitable material

Sample Weight 0.07 to 0.09 g

Analysis Time 7 minutes

Furnace Temperature 950°C

Flow Profile


#	Oxygen Flow	Time
1	Low	15
2	Low	15
3	Low	15
4	Medium	15
5	Medium	END

Procedure

1. Prepare the instrument by following the procedure as outlined in the operator's instruction manual (perform any required maintenance, perform leak checks, etc.).
2. Analyze gas blanks until instrument is stable. Reset the blank if necessary.
3. Switch to the above flow profile.
3. Analyze three to five 502-186 Tin Foil Cups containing ~0.07 to 0.09 g standard material and drift correct.
4. Weigh unknown sample into a 502-186 Tin Foil Cup, seal, load into carousel, and analyze.

Typical Results

Sample	% Carbon	% Hydrogen	% Nitrogen
Coke	88.94	0.96	1.74
	89.04	0.95	1.72
	88.70	0.94	1.76
	88.53	0.94	1.71
	89.16	0.95	1.71
	88.80	0.94	1.71
	89.07	0.94	1.71
	88.48	0.92	1.70
	88.45	0.94	1.71
	88.43	0.93	1.71
	88.43	0.93	1.72
Average	88.73	0.94	1.72
Std. Dev.	0.28	0.011	0.019



LECO Corporation • 3000 Lakeview Ave. • St. Joseph, MI 49085-2396
Phone: 800-292-6141 • Fax: 269-982-8977
info@leco.com • www.leco.com • ISO-9001 No. FM 24045

CHN-2000